

Table 2 Electrode locations and channel numbers of the layout in Fig. A.2. The column “cable” refers to cable pin number (Y: yellow, G: green).

Loc.	Channel	Cable	Loc.	Channel	Cable
FP1	1	Y1	AF7	33	G1
FP2	2	Y2	AF3	34	G2
F7	3	Y3	AFZ	35	G3
F3	4	Y4	AF4	36	G4
FZ	5	Y5	AF8	37	G5
F4	6	Y6	F5	38	G6
F8	7	Y7	F1	39	G7
T7	8	Y8	F2	40	G8
C3	9	Y9	F6	41	G9
CZ	10	Y10	FC3	42	G10
C4	11	Y11	FCZ	43	G11
T8	12	Y12	FC4	44	G12
P7	13	Y13	C5	45	G13
P3	14	Y14	C1	46	G14
PZ	15	Y15	C2	47	G15
P4	16	Y16	C6	48	G16
P8	17	Y17	CP3	49	G17
O1	18	Y18	CPZ	50	G18
O2	19	Y19	CP4	51	G19
FPZ	20	Y20	P5	52	G20
FC5	21	Y21	P1	53	G21
FC1	22	Y22	P2	54	G22
FC2	23	Y23	P6	55	G23
FC6	24	Y24	PO3	56	G24
CP5	25	Y25	PO4	57	G25
CP1	26	Y26	FT7	58	G26
CP2	27	Y27	FT8	59	G27
CP6	28	Y28	TP7	60	G28
PO5	29	Y29	TP8	61	G29
POZ	30	Y30	PO7	62	G30
PO6	31	Y31	PO8	63	G31
OZ	32	Y32	IZ	64	G32